



**Process for the fermentative preparation of L-amino acids
using strains of the Enterobacteriaceae family**

This invention relates to a process for the fermentative preparation of L-amino acids, in particular L-threonine, using strains of the Enterobacteriaceae family in which the poxB gene is attenuated.

Prior art

L-Amino acids, in particular L-threonine, are used in human medicine and in the pharmaceuticals industry, in the foodstuffs industry and very particularly in animal nutrition..(sic)

It is known to prepare L-amino acids by fermentation of strains of Enterobacteriaceae, in particular Escherichia coli (E. coli) and Serratia marcescens. Because of their great importance, work is constantly being undertaken to improve the preparation processes. Improvements to the process can relate to fermentation measures, such as e.g. stirring and supply of oxygen, or the composition of the nutrient media, such as e.g. the sugar concentration during the fermentation, or the working up to the product form, by e.g. ion exchange chromatography, or the intrinsic output properties of the microorganism itself.

Methods of mutagenesis, selection and mutant selection are used to improve the output properties of these microorganisms. Strains which are resistant to antimetabolites, such as e.g. the threonine analogue α -amino- β -hydroxyvaleric acid (AHV), or are auxotrophic for metabolites of regulatory importance and produce L-amino acid, such as e.g. L-threonine, are obtained in this manner.

Methods of the recombinant DNA technique have also been employed for some years for improving the strain of strains of the Enterobacteriaceae family which produce L-amino acids, by amplifying individual amino acid

biosynthesis genes and investigating the effect on the production.

Object of the invention

The inventors had the object of providing new measures for improved fermentative preparation of L-amino acids, in particular L-threonine.

Description of the invention

The invention provides a process for the fermentative preparation of L-amino acids, in particular L-threonine, using microorganisms of the Enterobacteriaceae family which, in particular, already produce L-threonine and in which the nucleotide sequence which codes for the enzyme pyruvate oxidase (EC 1.2.2.2) (poxB gene) is attenuated.

The term "attenuation" in this connection describes the reduction or elimination of the intracellular activity of one or more enzymes (proteins) in a microorganism which are coded by the corresponding DNA, for example by using a weak promoter or a gene or allele which codes for a corresponding enzyme with a low activity or inactivates the corresponding enzyme (protein) or gene, and optionally combining these measures.

The process comprises carrying out the following steps:

- a) fermentation of microorganisms of the Enterobacteriaceae family in which at least the poxB gene is attenuated,
- b) concentration of the corresponding L-amino acid in the medium or in the cells of the microorganisms of the Enterobacteriaceae family, and
- c) isolation of the desired L-amino acid.

The microorganisms which the present invention provides can produce L-amino acids from glucose, sucrose, lactose, fructose, maltose, molasses, optionally starch, optionally

cellulose or from glycerol and ethanol. They are representatives of the Enterobacteriaceae family chosen from the genera Escherichia, Erwinia, Providencia and Serratia. The genera Escherichia and Serratia are preferred. Of the genus Escherichia the species Escherichia coli and of the genus Serratia the species Serratia marcescens are to be mentioned in particular.

Suitable strains, which produce L-threonine in particular, of the genus Escherichia, in particular of the species Escherichia coli, are, for example

Escherichia coli TF427
Escherichia coli H4578
Escherichia coli KY10935
Escherichia coli VNIIGenetika MG442
Escherichia coli VNIIGenetika M1
Escherichia coli VNIIGenetika 472T23
Escherichia coli BKIIM B-3996
Escherichia coli kat 13
Escherichia coli KCCM-10132

Suitable L-threonine-producing strains of the genus Serratia, in particular of the species Serratia marcescens, are, for example

Serratia marcescens HNr21
Serratia marcescens TLr156
Serratia marcescens T2000

Strains from the Enterobacteriaceae family which produce L-threonine preferably have, inter alia, one or more genetic or phenotypic features chosen from the group consisting of: resistance to α -amino- β -hydroxyvaleric acid, resistance to thialysine, resistance to ethionine, resistance to α -methylserine, resistance to diaminosuccinic acid, resistance to α -aminobutyric acid, resistance to borrelidin, resistance to rifampicin, resistance to valine analogues, such as, for example, valine hydroxamate,

resistance to purine analogues, such as, for example, 6-dimethylaminopurine, a need for L-methionine, optionally a partial and compensatable need for L-isoleucine, a need for meso-diaminopimelic acid, auxotrophy in respect of

5 threonine-containing dipeptides, resistance to L-threonine, resistance to L-homoserine, resistance to L-lysine, resistance to L-methionine, resistance to L-glutamic acid, resistance to L-aspartate, resistance to L-leucine, resistance to L-phenylalanine, resistance to L-serine,

10 resistance to L-cysteine, resistance to L-valine, sensitivity to fluoropyruvate, defective threonine dehydrogenase, optionally a capacity for sucrose utilization, enhancement of the threonine operon, enhancement of homoserine dehydrogenase I-aspartate kinase

15 I, preferably of the feedback-resistant form, enhancement of homoserine kinase, enhancement of threonine synthase, enhancement of aspartate kinase, optionally of the feedback-resistant form, enhancement of aspartate semialdehyde dehydrogenase, enhancement of phosphoenol

20 pyruvate carboxylase, optionally of the feedback-resistant form, enhancement of phosphoenol pyruvate synthase, enhancement of transhydrogenase, enhancement of the RhtB gene product, enhancement of the RhtC gene product, enhancement of the YfiK gene product, enhancement of a

25 pyruvate carboxylase, and attenuation of acetic acid formation.

It has been found that microorganisms of the Enterobacteriaceae family produce L-amino acids, in particular L-threonine, in an improved manner after

30 attenuation, in particular-elimination, of the poxB gene, which codes for pyruvate oxidase (EC number 1.2.2.2).

It has furthermore been found that microorganisms of the Enterobacteriaceae family form lower concentrations of the undesirable by-product acetic acid after attenuation, in

particular elimination, of the poxB gene, which codes for pyruvate oxidase (EC number 1.2.2.2).

The nucleotide sequence of the poxB gene of Escherichia coli has been published by Grabau and Cronan (Nucleic Acids Research. 14 (13), 5449-5460 (1986)) and can also be found from the genome sequence of Escherichia coli published by Blattner et al. (Science 277, 1453 - 1462 (1997), under Accession Number AE000188. The nucleotide sequence of the poxB gene of Escherichia coli is shown in SEQ ID No. 1 and the amino acid sequence of the associated gene product is shown in SEQ ID No. 2.

The poxB genes described in the text references mentioned can be used according to the invention. Alleles of the poxB gene which result from the degeneracy of the genetic code or due to "sense mutations" of neutral function can furthermore be used.

To achieve an attenuation, for example, expression of the poxB gene or the catalytic properties of the enzyme protein can be reduced or eliminated. The two measures can optionally be combined.

The reduction in gene expression can take place by suitable culturing, by genetic modification (mutation) of the signal structures of gene expression or also by the antisense-RNA technique. Signal structures of gene expression are, for example, repressor genes, activator genes, operators, promoters, attenuators, ribosome binding sites, the start codon and terminators. The expert can find information in this respect, inter alia, for example, in Jensen and Hammer (Biotechnology and Bioengineering 58: 191-195 (1998)), in Carrier and Keasling (Biotechnology Progress 15, 58-64 (1999), Franch and Gerdes (Current Opinion in Microbiology 3, 159-164 (2000)) and in known textbooks of genetics and molecular biology, such as, for example, the textbook of Knippers ("Molekulare Genetik [Molecular Genetics]", 6th edition, Georg Thieme Verlag, Stuttgart, Germany, 1995) or

that of Winnacker ("Gene und Klone [Genes and Clones]", VCH Verlagsgesellschaft, Weinheim, Germany, 1990).

Mutations which lead to a change or reduction in the catalytic properties of enzyme proteins are known from the prior art. Examples which may be mentioned are the works of Qiu and Goodman (Journal of Biological Chemistry 272: 8611-8617 (1997)), Yano et al. (Proceedings of the National Academy of Sciences, USA 95, 5511-5515 (1998)), Wentz and Schachmann (Journal of Biological Chemistry 266, 20833-20839 (1991)). Summarizing descriptions can be found in known textbooks of genetics and molecular biology, such as e.g. that by Hagemann ("Allgemeine Genetik [General Genetics]", Gustav Fischer Verlag, Stuttgart, 1986).

Possible mutations are transitions, transversions, insertions and deletions. Depending on the effect of the amino acid exchange on the enzyme activity, "missense mutations" or "nonsense mutations" are referred to. Insertions or deletions of at least one base pair in a gene lead to "frame shift mutations", which lead to incorrect amino acids being incorporated or translation being interrupted prematurely. Deletions of several codons typically lead to a complete loss of the enzyme activity. Instructions on generation of such mutations are prior art and can be found in known textbooks of genetics and molecular biology, such as e.g. the textbook by Knippers ("Molekulare Genetik [Molecular Genetics]", 6th edition, Georg Thieme Verlag, Stuttgart, Germany, 1995), that by Winnacker ("Gene und Klone [Genes and Clones]", VCH Verlagsgesellschaft, Weinheim, Germany, 1990) or that by Hagemann ("Allgemeine Genetik [General Genetics]", Gustav Fischer Verlag, Stuttgart, 1986).

An example of a plasmid with the aid of which the poxB gene of Escherichia coli can be attenuated, in particular eliminated, by position-specific mutagenesis is the plasmid pMAK705 Δ poxB (figure 1). In addition to residues of polylinker sequences, it contains only a part of the 5' and

a part of the 3' region of the poxB gene. A 340 bp long section of the coding region is missing (deletion). The sequence of this DNA which can be employed for mutagenesis of the poxB gene is shown in SEQ ID No. 3.

- 5 The deletion mutation of the poxB gene can be incorporated into suitable strains by gene or allele replacement.

A conventional method is the method, described by Hamilton et al. (Journal of Bacteriology 174, 4617 - 4622 (1989)), of gene replacement with the aid of a conditionally
10 replicating pSC101 derivative pMAK705. Other methods described in the prior art, such as, for example, those of Martinez-Morales et al. (Journal of Bacteriology 1999, 7143-7148 (1999)) or those of Boyd et al. (Journal of Bacteriology 182, 842-847 (2000)), can likewise be used.

- 15 After replacement has taken place, the strain in question contains the form of the Δ poxB allele shown in SEQ ID No. 4, which is also provided by the invention.

It is also possible to transfer mutations in the poxB gene or mutations which affect expression of the poxB gene into
20 various strains by conjugation or transduction.

It may furthermore be advantageous for the production of L-amino acids, in particular L-threonine, with strains of the Enterobacteriaceae family to enhance one or more enzymes of the known threonine biosynthesis pathway or enzymes of
25 anaplerotic metabolism or enzymes for the production of reduced nicotinamide adenine dinucleotide phosphate, in addition to the attenuation of the poxB gene.

The term "enhancement" in this connection describes the increase in the intracellular activity of one or more
30 enzymes or proteins in a microorganism which are coded by the corresponding DNA, for example by increasing the number of copies of the gene or genes, using a potent promoter or a gene which codes for a corresponding enzyme or protein

with a high activity, and optionally combining these measures.

Thus, for example, one or more genes chosen from the group consisting of

- 5 • the thrABC operon which codes for aspartate kinase, homoserine dehydrogenase, homoserine kinase and threonine synthase (US-A-4,278,765),
 - the pyc gene which codes for pyruvate carboxylase (DE-A-19 831 609),
 - 10 • the pps gene which codes for phosphoenol pyruvate synthase (Molecular and General Genetics 231:332 (1992)),
 - the ppc gene which codes for phosphoenol pyruvate carboxylase (Gene 31:279-283 (1984)),
 - 15 • the pntA and pntB genes which code for transhydrogenase (European Journal of Biochemistry 158:647-653 (1986)),
 - the rhtB gene which imparts homoserine resistance (EP-A-0 994 190),
 - the mgo gene which codes for malate:quinone
20 oxidoreductase (DE 100 348 33.5),
 - the rhtC gene which imparts threonine resistance (EP-A-1 013 765), and
 - the thrE gene of Corynebacterium glutamicum which codes for threonine export (DE 100 264 94.8)
- 25 can be enhanced, in particular over-expressed, at the same time.

It may furthermore be advantageous for the production of L-amino acids, in particular threonine, in addition to the attenuation of the poxB gene, for one or more genes chosen
30 from the group consisting of

- the tdh gene which codes for threonine dehydrogenase (Ravnikar and Somerville, Journal of Bacteriology 169, 4716-4721 (1987)),
- the mdh gene which codes for malate dehydrogenase (E.C. 1.1.1.37) (Vogel et al., Archives in Microbiology 149, 36-42 (1987)),
- the gene product of the open reading frame (orf) yjfA (Accession Number AAC77180 of the National Center for Biotechnology Information (NCBI, Bethesda, MD, USA),
- the gene product of the open reading frame (orf) ytfP (Accession Number AAC77179 of the National Center for Biotechnology Information (NCBI, Bethesda, MD, USA) and
- the pckA gene which codes for the enzyme phosphoenolpyruvate carboxykinase (Medina et al. (Journal of Bacteriology 172, 7151-7156 (1990))

to be attenuated, in particular eliminated or reduced in expression.

In addition to attenuation of the poxB gene it may furthermore be advantageous for the production of L-amino acids, in particular L-threonine, to eliminate undesirable side reactions (Nakayama: "Breeding of Amino Acid Producing Microorganisms", in: Overproduction of Microbial Products, Krumphanzl, Sikyta, Vanek (eds.), Academic Press, London, UK, 1982).

The microorganisms produced according to the invention can be cultured in the batch process (batch culture), the fed batch (feed process) or the repeated fed batch process (feed process). A summary of known culture methods are [sic] described in the textbook by Chmiel (Bioprozesstechnik 1. Einführung in die Bioverfahrenstechnik [Bioprocess Technology 1. Introduction to Bioprocess Technology (Gustav Fischer Verlag, Stuttgart, 1991)) or in the textbook by Storhas (Bioreaktoren und

periphere Einrichtungen [Bioreactors and Peripheral Equipment] (Vieweg Verlag, Braunschweig/Wiesbaden, 1994)).

The culture medium to be used must meet the requirements of the particular strains in a suitable manner. Descriptions of culture media for various microorganisms are contained in the handbook "Manual of Methods for General Bacteriology" of the American Society for Bacteriology (Washington D.C., USA, 1981).

Sugars and carbohydrates, such as e.g. glucose, sucrose, lactose, fructose, maltose, molasses, starch and optionally cellulose, oils and fats, such as e.g. soya oil, sunflower oil, groundnut oil and coconut fat, fatty acids, such as e.g. palmitic acid, stearic acid and linoleic acid, alcohols, such as e.g. glycerol and ethanol, and organic acids, such as e.g. acetic acid, can be used as the source of carbon. These substances can be used individually or as a mixture.

Organic nitrogen-containing compounds, such as peptones, yeast extract, meat extract, malt extract, corn steep liquor, soya bean flour and urea, or inorganic compounds, such as ammonium sulfate, ammonium chloride, ammonium phosphate, ammonium carbonate and ammonium nitrate, can be used as the source of nitrogen. The sources of nitrogen can be used individually or as a mixture.

Phosphoric acid, potassium dihydrogen phosphate or dipotassium hydrogen phosphate or the corresponding sodium-containing salts can be used as the source of phosphorus. The culture medium must furthermore comprise salts of metals, such as e.g. magnesium sulfate or iron sulfate, which are necessary for growth. Finally, essential growth substances, such as amino acids and vitamins, can be employed in addition to the abovementioned substances. Suitable precursors can moreover be added to the culture medium. The starting substances mentioned can be added to the culture in the form of a single batch, or can be fed in during the culture in a suitable manner.

Basic compounds, such as sodium hydroxide, potassium hydroxide, ammonia or aqueous ammonia, or acid compounds, such as phosphoric acid or sulfuric acid, can be employed in a suitable manner to control the pH. Antifoams, such as
5 e.g. fatty acid polyglycol esters, can be employed to control the development of foam. Suitable substances having a selective action, e.g. antibiotics, can be added to the medium to maintain the stability of plasmids. To maintain aerobic conditions, oxygen or oxygen-containing
10 gas mixtures, such as e.g. air, are introduced into the culture. The temperature of the culture is usually 25°C to 45°C, and preferably 30°C to 40°C. Culturing is continued until a maximum of L-amino acids or L-threonine has formed. This target is usually reached within 10 hours to 160
15 hours.

The analysis of L-amino acids can be carried out by anion exchange chromatography with subsequent ninhydrin derivatization, as described by Spackman et al. (Analytical Chemistry, 30, (1958), 1190), or it can take place by
20 reversed phase HPLC as described by Lindroth et al. (Analytical Chemistry (1979) 51: 1167-1174).

A pure culture of the Escherichia coli K-12 strain DH5 α /pMAK705 was deposited as DSM 13720 on 12th September 2000 at the Deutsche Sammlung für Mikroorganismen und
25 Zellkulturen (DSMZ = German Collection of Microorganisms and Cell Cultures, Braunschweig, Germany) in accordance with the Budapest Treaty.

A pure culture of the Escherichia coli K-12 strain MG442 Δ poxB was deposited as DSM 13762 on 2nd October 2000
30 at the Deutsche Sammlung für Mikroorganismen und Zellkulturen (DSMZ = German Collection of Microorganisms and Cell Cultures, Braunschweig, Germany) in accordance with the Budapest Treaty.

The process according to the invention is used for the
35 fermentative preparation of L-amino acids, such as e.g.

L-threonine, L-isoleucine, L-valine, L-methionine, L-homoserine and L-lysine, in particular L-threonine.

The present invention is explained in more detail in the following with the aid of embodiment examples.

- 5 The isolation of plasmid DNA from *Escherichia coli* and all techniques of restriction, Klenow and alkaline phosphatase treatment are carried out by the method of Sambrook et al. (Molecular cloning - A laboratory manual (1989) Cold Spring Harbour Laboratory Press). Unless described otherwise, the
10 transformation of *Escherichia coli* is carried out by the method of Chung et al. (Proceedings of the National Academy of Sciences of the United States of America USA (1989) 86: 2172-2175).

The incubation temperature for the preparation of strains
15 and transformants is 37°C. Temperatures of 30°C and 44°C are used in the gene replacement method of Hamilton et. al.

Example 1

Construction of the deletion mutation of the *poxB* gene

Parts of the 5' and 3' region of the *poxB* gene are
20 amplified from *Escherichia coli* K12 using the polymerase chain reaction (PCR) and synthetic oligonucleotides. Starting from the nucleotide sequence of the *poxB* gene in *E. coli* K12 MG1655 (SEQ ID No. 1), the following PCR primers are synthesized (MWG Biotech, Ebersberg, Germany):

25 *poxB*'5'-1: 5' - CTGAACGGTCTTAGTGACAG - 3'

poxB'5'-2: 5' - AGGCCTGGAATAACGCAGCAGTTG - 3'

poxB'3'-1: 5' - CTGCGTGCATTGCTTCCATTG - 3'

poxB'3'-2: 5' - GCCAGTTCGATCACTTCATCAC - 3'

The chromosomal *E. coli* K12 MG1655 DNA employed for the PCR
30 is isolated according to the manufacturers instructions with "Qiagen Genomic-tips 100/G" (QIAGEN, Hilden, Germany).

A DNA fragment approx. 500 base pairs (bp) in size from the 5' region of the *poxB* gene (called *poxB1*) and a DNA fragment approx. 750 bp in size from the 3' region of the *poxB* gene (called *poxB2*) can be amplified with the specific primers under standard PCR conditions (Innis et al. (1990) PCR Protocols. A Guide to Methods and Applications, Academic Press) with Taq-DNA polymerase (Gibco-BRL, Eggenstein, Germany). The PCR products are each ligated with the vector pCR2.1TOPO (TOPO TA Cloning Kit, Invitrogen, Groningen, The Netherlands) in accordance with the manufacturers instructions and transformed into the *E. coli* strain TOP10F'.

Selection of plasmid-carrying cells takes place on LB agar, to which 50 µg/ml ampicillin are added. After isolation of the plasmid DNA, the vector pCR2.1TOPO*poxB1* is cleaved with the restriction enzymes *Ecl136II* and *XbaI* and, after separation in 0.8% agarose gel, the *poxB1* fragment is isolated with the aid of the QIAquick Gel Extraction Kit (QIAGEN, Hilden, Germany). After isolation of the plasmid DNA the vector pCR2.1TOPO*poxB2* is cleaved with the enzymes *EcoRV* and *XbaI* and ligated with the *poxB1* fragment isolated. The *E. coli* strain DH5α is transformed with the ligation batch and plasmid-carrying cells are selected on LB agar, to which 50 µg/ml ampicillin is added. After isolation of the plasmid DNA those plasmids in which the mutagenic DNA sequence shown in SEQ ID No. 3 is cloned are detected by control cleavage with the enzymes *HindIII* and *XbaI*. One of the plasmids is called pCR2.1TOPOΔ*poxB*.

Example 2

30 Construction of the replacement vector pMAK705Δ*poxB*

The *poxB* allele described in example 1 is isolated from the vector pCR2.1TOPOΔ*poxB* after restriction with the enzymes *HindIII* and *XbaI* and separation in 0.8% agarose gel, and ligated with the plasmid pMAK705 (Hamilton et al. (1989) Journal of Bacteriology 174, 4617 - 4622), which has been digested with the enzymes *HindIII* and *XbaI*. The ligation

batch is transformed in DH5 α and plasmid-carrying cells are selected on LB agar, to which 20 μ g/ml chloramphenicol is added. Successful cloning is demonstrated after isolation of the plasmid DNA and cleavage with the enzymes HindIII and XbaI. The replacement vector formed, pMAK705 Δ poxB (= pMAK705 Δ poxB), is shown in figure 1.

Example 3

Position-specific mutagenesis of the poxB gene in the E. coli strain MG442

10 The L-threonine-producing E. coli strain MG442 is described in the patent specification US-A- 4,278,765 and deposited as CMIM B-1628 at the Russian National Collection for Industrial Microorganisms (VKPM, Moscow, Russia).

For replacement of the chromosomal poxB gene with the
15 plasmid-coded deletion construct, MG442 is transformed with the plasmid pMAK705 Δ poxB. The gene replacement is carried out by the selection method described by Hamilton et al. (1989) Journal of Bacteriology 174, 4617 - 4622) and is verified by standard PCR methods (Innis et al. (1990) PCR
20 Protocols. A Guide to Methods and Applications, Academic Press) with the following oligonucleotide primers:

poxB'5'-1: 5' - CTGAACGGTCTTAGTGACAG - 3'

poxB'3'-2: 5' - GCCAGTTCGATCACTTCATCAC -3'

The strain obtained is called MG442 Δ poxB.

25 Example 4

Preparation of L-threonine with the strain MG442 Δ poxB

MG442 Δ poxB is multiplied on minimal medium with the following composition: 3.5 g/l Na₂HPO₄*2H₂O, 1.5 g/l KH₂PO₄, 1 g/l NH₄Cl, 0.1 g/l MgSO₄*7H₂O, 2 g/l glucose, 20 g/l agar.
30 The formation of L-threonine is checked in batch cultures of 10 ml contained in 100 ml conical flasks. For this, 10 ml of preculture medium of the following composition:

2 g/l yeast extract, 10 g/l $(\text{NH}_4)_2\text{SO}_4$, 1 g/l KH_2PO_4 , 0.5 g/l $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$, 15 g/l CaCO_3 , 20 g/l glucose are inoculated and the batch is incubated for 16 hours at 37°C and 180 rpm on an ESR incubator from Kühner AG (Birsfelden, Switzerland).
5 250 µl of this preculture are transinoculated into 10 ml of production medium (25 g/l $(\text{NH}_4)_2\text{SO}_4$, 2 g/l KH_2PO_4 , 1 g/l $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$, 0.03 g/l $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$, 0.018 g/l $\text{MnSO}_4 \cdot \text{H}_2\text{O}$, 30 g/l CaCO_3 , 20 g/l glucose) and the batch is incubated for 48 hours at 37°C. After the incubation the optical density
10 (OD) of the culture suspension is determined with an LP2W photometer from Dr. Lange (Berlin, Germany) at a measurement wavelength of 660 nm.

The concentration of L-threonine formed is then determined in the sterile-filtered culture supernatant with an amino
15 acid analyzer from Eppendorf-BioTronik (Hamburg, Germany) by ion exchange chromatography and post-column reaction with ninhydrin detection.

The result of the experiment is shown in table 1.

Table 1

Strain	OD (660 nm)	L-Threonine g/l
MG442	6.0	1.5
MG442 Δ poxB	4.9	2.6

The following figure is attached:

- Figure 1: pMAK705 Δ poxB (= pMAK705deltapoxB)

5 The length data are to be understood as approx. data. The abbreviations and designations used have the following meaning:

- cat: chloramphenicol resistance gene
- rep-ts: temperature-sensitive replication region of
10 the plasmid pSC101
- poxB1: part of the 5' region of the poxB gene
- poxB2: part of the 3' region of the poxB gene

The abbreviations for the restriction enzymes have the following meaning

- 15 • BamHI: restriction endonuclease from *Bacillus amyloliquefaciens*
- BglII: restriction endonuclease from *Bacillus globigii*
- ClaI: restriction endonuclease from *Caryophanon latum*
- Ecl136II restriction endonuclease from *Enterobacter*
20 *cloacae* RFL136 (= Ecl136)
- EcoRI: restriction endonuclease from *Escherichia coli*
- EcoRV: restriction endonuclease from *Escherichia coli*

- HindIII: restriction endonuclease from *Haemophilus influenzae*
- KpnI: restriction endonuclease from *Klebsiella pneumoniae*
- 5 • PstI: restriction endonuclease from *Providencia stuartii*
- PvuI: restriction endonuclease from *Proteus vulgaris*
- SmaI: restriction endonuclease from *Serratia marcescens*
- 10 • XbaI: restriction endonuclease from *Xanthomonas badrii*
- XhoI: restriction endonuclease from *Xanthomonas holcicola*

SEQUENCE PROTOCOL

<110> Degussa-Hüls AG

5 <120> Process for the fermentative preparation of
L-amino acids using strains of the
Enterobacteriaceae family.

10 <130> 000613 BT

<140>
<141>

15 <160> 4

<170> PatentIn Ver. 2.1

20 <210> 1
<211> 1719
<212> DNA
<213> Escherichia coli

25 <220>
<221> CDS
<222> (1)..(1716)
<223> poxB

<400> 1

30	atg	aaa	caa	acg	gtt	gca	gct	tat	atc	gcc	aaa	aca	ctc	gaa	tcg	gca	48
	Met	Lys	Gln	Thr	Val	Ala	Ala	Tyr	Ile	Ala	Lys	Thr	Leu	Glu	Ser	Ala	
	1				5					10					15		
35	ggg	gtg	aaa	cgc	atc	tgg	gga	gtc	aca	ggc	gac	tct	ctg	aac	ggt	ctt	96
	Gly	Val	Lys	Arg	Ile	Trp	Gly	Val	Thr	Gly	Asp	Ser	Leu	Asn	Gly	Leu	
				20				25						30			
40	agt	gac	agt	ctt	aat	cgc	atg	ggc	acc	atc	gag	tgg	atg	tcc	acc	cgc	144
	Ser	Asp	Ser	Leu	Asn	Arg	Met	Gly	Thr	Ile	Glu	Trp	Met	Ser	Thr	Arg	
				35				40					45				
45	cac	gaa	gaa	gtg	gcg	gcc	ttt	gcc	gct	ggc	gct	gaa	gca	caa	ctt	agc	192
	His	Glu	Glu	Val	Ala	Ala	Phe	Ala	Ala	Gly	Ala	Glu	Ala	Gln	Leu	Ser	
				50			55					60					
50	gga	gaa	ctg	gcg	gtc	tgc	gcc	gga	tcg	tgc	ggc	ccc	ggc	aac	ctg	cac	240
	Gly	Glu	Leu	Ala	Val	Cys	Ala	Gly	Ser	Cys	Gly	Pro	Gly	Asn	Leu	His	
						70					75					80	
55	tta	atc	aac	ggc	ctg	ttc	gat	tgc	cac	cgc	aat	cac	gtt	ccg	gta	ctg	288
	Leu	Ile	Asn	Gly	Leu	Phe	Asp	Cys	His	Arg	Asn	His	Val	Pro	Val	Leu	
					85					90					95		
60	gcg	att	gcc	gct	cat	att	ccc	tcc	agc	gaa	att	ggc	agc	ggc	tat	ttc	336
	Ala	Ile	Ala	Ala	His	Ile	Pro	Ser	Ser	Glu	Ile	Gly	Ser	Gly	Tyr	Phe	
				100				105						110			
65	cag	gaa	acc	cac	cca	caa	gag	cta	ttc	cgc	gaa	tgt	agt	cac	tat	tgc	384
	Gln	Glu	Thr	His	Pro	Gln	Glu	Leu	Phe	Arg	Glu	Cys	Ser	His	Tyr	Cys	
				115				120					125				
70	gag	ctg	gtt	tcc	agc	ccg	gag	cag	atc	cca	caa	gta	ctg	gcg	att	gcc	432
	Glu	Leu	Val	Ser	Ser	Pro	Glu	Gln	Ile	Pro	Gln	Val	Leu	Ala	Ile	Ala	
				130			135					140					
75	atg	cgc	aaa	gcg	gtg	ctt	aac	cgt	ggc	gtt	tcg	gtt	gtc	gtg	tta	cca	480
	Met	Arg	Lys	Ala	Val	Leu	Asn	Arg	Gly	Val	Ser	Val	Val	Val	Leu	Pro	

	145		150		155		160										
5	ggc Gly	gac Asp	gtg Val	gcg Ala	tta Leu 165	aaa Lys	cct Pro	gcg Ala	cca Pro	gaa Glu 170	ggg Gly	gca Ala	acc Thr	atg Met	cac His 175	tgg Trp	528
10	tat Tyr	cat His	gcg Ala	cca Pro 180	caa Gln	cca Pro	gtc Val	gtg Val	acg Thr 185	ccg Pro	gaa Glu	gaa Glu	gaa Glu	gag Glu 190	tta Leu	cgc Arg	576
	aaa Lys	ctg Leu	gcg Ala 195	caa Gln	ctg Leu	ctg Leu	cgt Arg	tat Tyr 200	tcc Ser	agc Ser	aat Asn	atc Ile	gcc Ala 205	ctg Leu	atg Met	tgt Cys	624
15	ggc Gly	agc Ser 210	ggc Gly	tgc Cys	gcg Ala	ggg Gly	gcg Ala 215	cat His	aaa Lys	gag Glu	tta Leu	gtt Val 220	gag Glu	ttt Phe	gcc Ala	ggg Gly	672
20	aaa Lys 225	att Ile	aaa Lys	gcg Ala	cct Pro	att Ile 230	gtt Val	cat His	gcc Ala	ctg Leu	cgc Arg 235	ggt Gly	aaa Lys	gaa Glu	cat His	gtc Val 240	720
25	gaa Glu	tac Tyr	gat Asp	aat Asn 245	ccg Pro	tat Tyr	gat Asp	gtt Val	gga Gly	atg Met 250	acc Thr	ggg Gly	tta Leu	atc Ile	ggc Gly 255	ttc Phe	768
30	tcg Ser	tca Ser	ggt Gly	ttc Phe 260	cat His	acc Thr	atg Met	atg Met	aac Asn 265	gcc Ala	gac Asp	acg Thr	tta Leu	gtg Val 270	cta Leu	ctc Leu	816
	ggc Gly	acg Thr	caa Gln 275	ttt Phe	ccc Pro	tac Tyr	cgc Arg	gcc Ala 280	ttc Phe	tac Tyr	ccg Pro	acc Thr	gat Asp 285	gcc Ala	aaa Lys	atc Ile	864
35	att Ile	cag Gln 290	att Ile	gat Asp	atc Ile	aac Asn	cca Pro 295	gcc Ala	agc Ser	atc Ile	ggc Gly	gct Ala 300	cac His	agc Ser	aag Lys	gtg Val	912
40	gat Asp 305	atg Met	gca Ala	ctg Leu	gtc Val	ggc Gly 310	gat Asp	atc Ile	aag Lys	tcg Ser	act Thr 315	ctg Leu	cgt Arg	gca Ala	ttg Leu	ctt Leu 320	960
45	cca Pro	ttg Leu	gtg Val	gaa Glu	gaa Glu 325	aaa Lys	gcc Ala	gat Asp	cgc Arg	aag Lys 330	ttt Phe	ctg Leu	gat Asp	aaa Lys	gcg Ala 335	ctg Leu	1008
50	gaa Glu	gat Asp	tac Tyr	cgc Arg 340	gac Asp	gcc Ala	cgc Arg	aaa Lys	ggg Gly 345	ctg Leu	gac Asp	gat Asp	tta Leu	gct Ala 350	aaa Lys	ccg Pro	1056
	agc Ser	gag Glu	aaa Lys 355	gcc Ala	att Ile	cac His	ccg Pro	caa Gln 360	tat Tyr	ctg Leu	gcg Ala	cag Gln 365	caa Gln	att Ile	agt Ser	cat His	1104
55	ttt Phe	gcc Ala 370	gcc Ala	gat Asp	gac Asp	gct Ala	att Ile 375	ttc Phe	acc Thr	tgt Cys	gac Asp	gtt Val 380	ggt Gly	acg Thr	cca Pro	acg Thr	1152
60	gtg Val 385	tgg Trp	gcg Ala	gca Ala	cgt Arg	tat Tyr 390	cta Leu	aaa Lys	atg Met	aac Asn	ggc Gly 395	aag Lys	cgt Arg	cgc Arg	ctg Leu	tta Leu 400	1200
65	ggt Gly	tcg Ser	ttt Phe	aac Asn 405	cac His	ggt Gly	tcg Ser	atg Met	gct Ala	aac Asn 410	gcc Ala	atg Met	ccg Pro	cag Gln 415	gcg Ala	ctg Leu	1248

	ggt gcg cag gcg aca gag cca gaa cgt cag gtg gtc gcc atg tgc ggc	1296
	Gly Ala Gln Ala Thr Glu Pro Glu Arg Gln Val Val Ala Met Cys Gly	
	420 425 430	
5	gat ggc ggt ttt agc atg ttg atg ggc gat ttc ctc tca gta gtg cag	1344
	Asp Gly Gly Phe Ser Met Leu Met Gly Asp Phe Leu Ser Val Val Gln	
	435 440 445	
10	atg aaa ctg cca gtg aaa att gtc gtc ttt aac aac agc gtg ctg ggc	1392
	Met Lys Leu Pro Val Lys Ile Val Val Phe Asn Asn Ser Val Leu Gly	
	450 455 460	
15	ttt gtg gcg atg gag atg aaa gct ggt ggc tat ttg act gac ggc acc	1440
	Phe Val Ala Met Glu Met Lys Ala Gly Gly Tyr Leu Thr Asp Gly Thr	
	465 470 475 480	
20	gaa cta cac gac aca aac ttt gcc cgc att gcc gaa gcg tgc ggc att	1488
	Glu Leu His Asp Thr Asn Phe Ala Arg Ile Ala Glu Ala Cys Gly Ile	
	485 490 495	
	acg ggt atc cgt gta gaa aaa gcg tct gaa gtt gat gaa gcc ctg caa	1536
	Thr Gly Ile Arg Val Glu Lys Ala Ser Glu Val Asp Glu Ala Leu Gln	
	500 505 510	
25	cgc gcc ttc tcc atc gac ggt ccg gtg ttg gtg gat gtg gtg gtc gcc	1584
	Arg Ala Phe Ser Ile Asp Gly Pro Val Leu Val Asp Val Val Val Ala	
	515 520 525	
30	aaa gaa gag tta gcc att cca ccg cag atc aaa ctc gaa cag gcc aaa	1632
	Lys Glu Glu Leu Ala Ile Pro Pro Gln Ile Lys Leu Glu Gln Ala Lys	
	530 535 540	
35	ggt ttc agc ctg tat atg ctg cgc gca atc atc agc gga cgc ggt gat	1680
	Gly Phe Ser Leu Tyr Met Leu Arg Ala Ile Ile Ser Gly Arg Gly Asp	
	545 550 555 560	
40	gaa gtg atc gaa ctg gcg aaa aca aac tgg cta agg taa	1719
	Glu Val Ile Glu Leu Ala Lys Thr Asn Trp Leu Arg	
	565 570	
	<210> 2	
	<211> 572	
	<212> PRT	
45	<213> Escherichia coli	
	<400> 2	
50	Met Lys Gln Thr Val Ala Ala Tyr Ile Ala Lys Thr Leu Glu Ser Ala	
	1 5 10 15	
	Gly Val Lys Arg Ile Trp Gly Val Thr Gly Asp Ser Leu Asn Gly Leu	
	20 25 30	
55	Ser Asp Ser Leu Asn Arg Met Gly Thr Ile Glu Trp Met Ser Thr Arg	
	35 40 45	
	His Glu Glu Val Ala Ala Phe Ala Ala Gly Ala Glu Ala Gln Leu Ser	
	50 55 60	
60		

	Gly	Glu	Leu	Ala	Val	Cys	Ala	Gly	Ser	Cys	Gly	Pro	Gly	Asn	Leu	His	
	65					70					75					80	
5	Leu	Ile	Asn	Gly	Leu	Phe	Asp	Cys	His	Arg	Asn	His	Val	Pro	Val	Leu	
					85					90					95		
	Ala	Ile	Ala	Ala	His	Ile	Pro	Ser	Ser	Glu	Ile	Gly	Ser	Gly	Tyr	Phe	
				100					105					110			
10	Gln	Glu	Thr	His	Pro	Gln	Glu	Leu	Phe	Arg	Glu	Cys	Ser	His	Tyr	Cys	
			115					120					125				
	Glu	Leu	Val	Ser	Ser	Pro	Glu	Gln	Ile	Pro	Gln	Val	Leu	Ala	Ile	Ala	
15		130					135					140					
	Met	Arg	Lys	Ala	Val	Leu	Asn	Arg	Gly	Val	Ser	Val	Val	Val	Leu	Pro	
	145					150					155				160		
20	Gly	Asp	Val	Ala	Leu	Lys	Pro	Ala	Pro	Glu	Gly	Ala	Thr	Met	His	Trp	
					165					170					175		
	Tyr	His	Ala	Pro	Gln	Pro	Val	Val	Thr	Pro	Glu	Glu	Glu	Glu	Leu	Arg	
				180					185					190			
25	Lys	Leu	Ala	Gln	Leu	Leu	Arg	Tyr	Ser	Ser	Asn	Ile	Ala	Leu	Met	Cys	
			195					200					205				
	Gly	Ser	Gly	Cys	Ala	Gly	Ala	His	Lys	Glu	Leu	Val	Glu	Phe	Ala	Gly	
30		210					215					220					
	Lys	Ile	Lys	Ala	Pro	Ile	Val	His	Ala	Leu	Arg	Gly	Lys	Glu	His	Val	
	225					230					235				240		
35	Glu	Tyr	Asp	Asn	Pro	Tyr	Asp	Val	Gly	Met	Thr	Gly	Leu	Ile	Gly	Phe	
				245						250					255		
	Ser	Ser	Gly	Phe	His	Thr	Met	Met	Asn	Ala	Asp	Thr	Leu	Val	Leu	Leu	
				260					265					270			
40	Gly	Thr	Gln	Phe	Pro	Tyr	Arg	Ala	Phe	Tyr	Pro	Thr	Asp	Ala	Lys	Ile	
			275					280					285				
	Ile	Gln	Ile	Asp	Ile	Asn	Pro	Ala	Ser	Ile	Gly	Ala	His	Ser	Lys	Val	
45		290				295					300						
	Asp	Met	Ala	Leu	Val	Gly	Asp	Ile	Lys	Ser	Thr	Leu	Arg	Ala	Leu	Leu	
	305					310					315				320		
50	Pro	Leu	Val	Glu	Glu	Lys	Ala	Asp	Arg	Lys	Phe	Leu	Asp	Lys	Ala	Leu	
				325						330					335		
	Glu	Asp	Tyr	Arg	Asp	Ala	Arg	Lys	Gly	Leu	Asp	Asp	Leu	Ala	Lys	Pro	
				340					345					350			
55	Ser	Glu	Lys	Ala	Ile	His	Pro	Gln	Tyr	Leu	Ala	Gln	Gln	Ile	Ser	His	
			355					360					365				
	Phe	Ala	Ala	Asp	Asp	Ala	Ile	Phe	Thr	Cys	Asp	Val	Gly	Thr	Pro	Thr	
60		370					375					380					
	Val	Trp	Ala	Ala	Arg	Tyr	Leu	Lys	Met	Asn	Gly	Lys	Arg	Arg	Leu	Leu	
	385					390					395				400		

Gly Ser Phe Asn His Gly Ser Met Ala Asn Ala Met Pro Gln Ala Leu
 405 410 415
 5 Gly Ala Gln Ala Thr Glu Pro Glu Arg Gln Val Val Ala Met Cys Gly
 420 425 430
 Asp Gly Gly Phe Ser Met Leu Met Gly Asp Phe Leu Ser Val Val Gln
 435 440 445
 10 Met Lys Leu Pro Val Lys Ile Val Val Phe Asn Asn Ser Val Leu Gly
 450 455 460
 Phe Val Ala Met Glu Met Lys Ala Gly Gly Tyr Leu Thr Asp Gly Thr
 465 470 475 480
 15 Glu Leu His Asp Thr Asn Phe Ala Arg Ile Ala Glu Ala Cys Gly Ile
 485 490 495
 Thr Gly Ile Arg Val Glu Lys Ala Ser Glu Val Asp Glu Ala Leu Gln
 500 505 510
 Arg Ala Phe Ser Ile Asp Gly Pro Val Leu Val Asp Val Val Val Ala
 515 520 525
 25 Lys Glu Glu Leu Ala Ile Pro Pro Gln Ile Lys Leu Glu Gln Ala Lys
 530 535 540
 Gly Phe Ser Leu Tyr Met Leu Arg Ala Ile Ile Ser Gly Arg Gly Asp
 545 550 555 560
 30 Glu Val Ile Glu Leu Ala Lys Thr Asn Trp Leu Arg
 565 570
 35
 <210> 3
 <211> 1454
 <212> DNA
 <213> Escherichia coli
 40
 <220>
 <221> misc_feature
 <222> (1)..(1454)
 <223> Mutagenic DNA
 45
 <220>
 <221> misc_feature
 <222> (1)..(56)
 <223> Technical DNA/residues of the polylinker sequence
 50
 <220>
 <221> misc_feature
 <222> (57)..(577)
 <223> Part of the 5' region (poxB1) of the poxB gene
 55
 <220>
 <221> misc_feature
 <222> (578)..(646)
 <223> Technical DNA/residues of the polylinker sequence
 60
 <220>
 <221> misc_feature
 <222> (647)..(1398)
 <223> Part of the 3' region (poxB2) of the poxB gene
 65
 <220>

<221> misc_feature
 <222> (1399)..(1454)
 <223> Technical DNA/residues of the polylinker sequence

```

5  <400> 3
   ctagatgcat gctcgagcgg ccgccagtgt gatggatatc tgcagaattc gcccttctga 60
   acggtcttag tgacagtctt aatcgcatgg gcaccatcga gtggatgtcc acccgccacg 120
   aagaagtggc ggcctttgcc gctggcgctg aagcacaact tagcggagaa ctggcggtct 180
   gcgcccggatc gtgcggcccc ggcaacctgc acttaatcaa cggcctgttc gattgccacc 240
10  gcaatcacgt tccggtactg gcgattgccg ctcatattcc ctccagcgaa attggcagcg 300
   gctattttcca ggaaaccac ccacaagagc tattccgcga atgtagtcac tattgcgagc 360
   tggtttccag cccggagcag atcccacaag tactggcgat tgccatgcgc aaagcgggtg 420
   ttaaccgtgg cgtttcggtt gtcgtgttac caggcgacgt ggcgttaaaa cctgcgccag 480
   aaggggcaac catgcactgg tatcatgcgc cacaaccagt cgtgacgccg gaagaagaag 540
15  agttacgcaa actggcgcaa ctgctgcgtt attccaggcc taagggcgaa ttccagcaca 600
   ctggcggccg ttactagtgg atccgagatc tgcagaattc gcccttctgc gtgcattgct 660
   tccattggtg gaagaaaaag ccgatcgcaa gtttctggat aaagcgctgg aagattaccg 720
   cgacgcccgc aaagggctgg acgatttagc taaaccgagc gagaaagcca ttcacccgca 780
   atatctggcg cagcaaatta gtcattttgc cgccgatgac gctattttca cctgtgacgt 840
20  tggtagccca acggtgtggg cggcacgtta tctaaaaatg aacggcaagc gtcgcctgtt 900
   aggttcgttt aaccacggtt cgatggctaa cgccatgccg caggcgctgg gtgcgcaggc 960
   gacagagcca gaacgtcagg tggtcgccat gtgcggcgat ggcggtttta gcatgttgat 1020
   gggcgatttc ctctcagtag tgcagatgaa actgccagtg aaaattgtcg tctttaacaa 1080
   cagcgtgctg ggctttgtgg cgatggagat gaaagctggt ggctatttga ctgacggcac 1140
25  cgaactacac gacacaaact ttgcccgcac tgccgaagcg tgcggcatta cgggtatccg 1200
   tgtagaaaaa gcgtctgaag ttgatgaagc cctgcaacgc gccttctcca tcgacgggtc 1260
   ggtgttggtg gatgtggtgg tcgccaaaga agagttagcc attccaccgc agatcaaact 1320
   cgaacaggcc aaaggtttca gcctgtatat gctgcgcgca atcatcagcg gacgcggtga 1380
   tgaagtgatc gaactggcaa gggcgaattc cagcacactg gcggccgtta ctagtggatc 1440
30  cgagctcggt acca                                     1454

```

<210> 4
 <211> 1448
 35 <212> DNA
 <213> Escherichia coli

<220>
 <221> misc_feature
 40 <222> (1)..(3)
 <223> Start codon of the delta poxB allele

<220>
 <221> misc_feature
 45 <222> (1)..(605)
 <223> 5' region of the delta poxB allele

<220>
 <221> misc_feature
 50 <222> (606)..(674)
 <223> Technical DNA/residues of the polylinker sequence

<220>
 <221> misc_feature
 55 <222> (675)..(1445)
 <223> 3' region of the delta poxB allele

<220>
 <221> misc_feature
 60 <222> (1446)..(1448)
 <223> Stop codon of the delta poxB allele

```

   <400> 4
65  atgaaacaaa cggttgcagc ttatatcgcc aaaacactcg aatcggcagg ggtgaaacgc 60
   atctggggag tcacaggcga ctctctgaac ggtcttagtg acagtcttaa tcgcatgggc 120
   accatcgagt ggatgtccac ccgccacgaa gaagtggcgg cctttgccgc tggcgctgaa 180

```

	gcacaactta	gcggagaact	ggcgggtctgc	gccggatcgt	gcggccccgg	caacctgcac	240
	ttaatcaacg	gcctgttcga	ttgccaccgc	aatcacgttc	cggtagctggc	gattgccgct	300
	catattccct	ccagcgaaat	tggcagcggc	tatttccagg	aaacccaccc	acaagagcta	360
	ttccgcgaat	gtagtcacta	ttgcgagctg	gtttccagcc	cggagcagat	cccacaagta	420
5	ctggcgattg	ccatgcgcaa	agcgggtgctt	aaccgtggcg	tttcggttgt	cgtgttacca	480
	ggcgacgtgg	cgttaaaacc	tgcgccagaa	ggggcaacca	tgcactggta	tcatgcgcca	540
	caaccagtcg	tgacgccgga	agaagaagag	ttacgcaaac	tggcgcaact	gctgcgttat	600
	tccaggccta	agggcgaaat	ccagcacact	ggcggccggt	actagtggat	ccgagatctg	660
	cagaattcgc	ccttctgcgt	gcattgcttc	cattgggtgga	agaaaaagcc	gatcgcaagt	720
10	ttctggataa	agcgctggaa	gattaccgcg	acgcccgcaa	agggtctggac	gatttagcta	780
	aaccgagcga	gaaagccatt	cacccgcaat	atctggcgca	gcaaattagt	cattttgccg	840
	ccgatgacgc	tattttcacc	tgtgacgttg	gtacgccaac	ggtgtgggcg	gcacgttatc	900
	taaaaatgaa	cggcaagcgt	cgcctgttag	gttcgtttta	ccacggttcg	atggctaacg	960
	ccatgccgca	ggcgctgggt	gcgcaggcga	cagagccaga	acgtcagggtg	gtcgccatgt	1020
15	gcggcgatgg	cggtttttagc	atgttgatgg	gcgatttcct	ctcagtagtg	cagatgaaac	1080
	tgccagtgaa	aattgtcgtc	tttaacaaca	gcgtgctggg	ctttgtggcg	atggagatga	1140
	aagctggtgg	ctatttgact	gacggcaccg	aactacacga	cacaaacttt	gcccgcattg	1200
	ccgaagcgtg	cggcattacg	ggtatccgtg	tagaaaaagc	gtctgaagtt	gatgaagccc	1260
	tgcaacgcgc	cttctccatc	gacgggccgg	tgttggtgga	tgtggtggtc	gccaaagaag	1320
20	agttagccat	tccaccgcag	atcaaactcg	aacaggccaa	aggtttcagc	ctgtatatgc	1380
	tgcgcgcaat	catcagcgga	cgcggtgatg	aagtgatcga	actggcgaaa	acaaactggc	1440
	taaggtaa						1448